

## How does AcceleDent work?

The appliance, used daily for 20 minutes, applies patented technology to gently vibrate the teeth and surrounding bone. In a controlled clinical trial, the vibrations have proven to accelerate tooth movement by 38% - 50%.

## Why are microvibrations so important?

The rate of tooth movement greatly depends on the speed at which the bone surrounding the teeth changes. Studies have shown that orthodontic force coupled with microvibrations, generated by the AcceleDent System, safely accelerate tooth movement.

## Is AcceleDent safe?

AcceleDent has a good safety profile, as indicated by university studies in humans. AcceleDent is considered safe under normal conditions of use with a conventional fixed orthodontic appliance.

## If it accelerates tooth movement, will it cause more discomfort?

No. AcceleDent is well tolerated and does not increase the level of discomfort associated with orthodontic treatment.

*"It is based on sound biological and clinical research and on the same principles used in the most current treatment of osteoporosis."*

*Dr. Dubravko Pavlin - Orthodontic Specialist*

*"One of the most exciting things about the technology is that it can be used together with braces, regardless of the type of bite problems."*

*Dr. Ravi Anthony - Orthodontic Specialist*

## Indications for Use:

AcceleDent is an orthodontic accessory intended for use during orthodontic treatment. It is used in conjunction with orthodontic appliances such as braces and helps facilitate minor anterior tooth movement.

Ask your Doctor about  
AcceleDent and visit  
[www.acceleddent.com](http://www.acceleddent.com)

Orthodontist Contact information:



8275 El Rio St. - Suite 100  
Houston, TX 77054 USA

L02100  
405 - 0406 - 001A

November 2011

North America:  
+1 832.631.1660  
Outside North America:  
+44 (0) 20 3318 1915  
[www.orthoaccel.com](http://www.orthoaccel.com)  
Rx only  
For orthodontic applications

# THE *FAST* TRACK TO STRAIGHT TEETH



**AcceleDent**<sup>™</sup>  
System



*"It's easy to use and comfortable."*

Alyssa U. - 14 year old patient  
San Antonio, TX

The AcceleDent™ System is designed to **INCREASE THE RATE OF TOOTH MOVEMENT** in patients undergoing treatment with a conventional fixed orthodontic appliance, which may reduce treatment time.

The use of AcceleDent with other types of orthodontic appliances has not been studied.

*"The best point is it's easy to use. I am extremely pleased with my orthodontic treatment and feel that AcceleDent has helped me achieve the cosmetic result I hoped for."*

Nicole S. - Adult patient, London, UK

### Patient Benefits

- Accelerated tooth movement without compromised aesthetics
- Compatible with braces for both adults and adolescents
- Lightweight for comfort
- Hands-free design allows you to go about your daily activities during the 20 minute treatment

### AcceleDent Components

#### 1 Activator

The Activator generates the patented microvibration transmitted to the teeth. Featuring a hands-free design, the Activator can be used almost anywhere.

#### 2 Mouthpiece

The Mouthpiece is accurately sized by your orthodontist to provide you balanced fit and comfort. It effortlessly unsnaps from the Activator for convenient cleaning.

#### 3 Charging Port

Docking the Activator with the Charging Port uploads compliance information to the "Smart Display" allowing you and your orthodontist to monitor the use.

#### 4 Travel Shell

The Travel Shell is an easy and effective way to store and transport your AcceleDent.



*"This technology will revolutionize orthodontic treatment for our patients."*

Dr. Lawrence J. Levens - Orthodontic Specialist